- Ⓔ *Venturi/Ejector Suction Module*
- Ⓓ *Venturi/Ejektor-Absaugmodul*
- Ⓕ *Module d'aspiration éjecteur/venturi*
- Ⓔ *Módulo de succión con eyector/válvula venturi*
- Ⓘ *Modulo di aspirazione venturi/eiettore*
- Ⓔ *Venturi/Ejektorsugmodul*
- Ⓔ *Модуль всасывания /эжектор Вентури*
- Ⓔ *De venturi-/ejector vacuümmodule*

Ⓔ	<i>Instructions for Use</i>	2
Ⓓ	<i>Bedienungsanleitung</i>	8
Ⓕ	<i>Mode d'emploi</i>	14
Ⓔ	<i>Instrucciones de uso</i>	20
Ⓘ	<i>Istruzioni per l'uso</i>	26
Ⓔ	<i>Bruksanvisning</i>	32
Ⓔ	<i>Инструкция по эксплуатации</i>	38
Ⓔ	<i>Gebruiksaanwijzing</i>	44


Table of Contents

Safety Instructions	2
Intended use	2
Receiving/Inspection	2
User Responsibility	2
Definitions	3
Venturi/ejector suction module identification	4
Equipment Set Up	4
Pre-use Checkout Procedure	5
Cleaning	6
Troubleshooting	6
Warranty	7

Safety Instructions

This manual provides you with important information about the Venturi/Ejector Suction Module. To ensure the safe and proper use of this device, READ and UNDERSTAND all of the safety and operating instructions. IF YOU DO NOT UNDERSTAND THESE INSTRUCTIONS, OR HAVE ANY QUESTIONS, CONTACT YOUR SUPERVISOR, DEALER OR THE MANUFACTURER BEFORE ATTEMPTING TO USE THE DEVICE.

Intended Use

The venturi/ ejector suction module is an accessory designed to generate a vacuum supply through the use of pressurized medical air or oxygen. It is not a stand alone medical device. It should always be used with a vacuum regulator. The device is to be used with Ohio Medical vacuum regulators only.

Receiving / Inspection:

Remove product from package and inspect for damage. If product is damaged, DO NOT USE and contact your dealer or equipment provider.

User Responsibility

This product performs as explained in this manual as long as the assembly, use, repair and maintenance are properly followed according to our instructions. Periodic review of this device is recommended. If any damage or defects are present, the product should not be used. This includes parts that may have been altered, contaminated, worn or missing. If any of the above are noted, immediate repair / replacement is required. In compliance with the Ohio Medical Warranty, repair of this device is not to be performed by anyone other than a qualified professional and done in strict accordance to the written instructions provided by Ohio Medical. If this device is subject to improper maintenance, repair, use and/or abuse leading to malfunction of the device, replacement is the sole responsibility of the user.

123 A This alpha character indicates the year of product manufacture and when the serial number was assigned; "Y" = 1995, "Z" = 1996, "A" = 1997, etc. "I" and "O" are not used.

WARNINGS  **This device is to be used only by persons who have been adequately instructed in its use.**

- △ **Never use any petroleum based lubricants in an Oxygen environment, as these materials are highly combustible in the presence of Oxygen. The only Oxygen service lubricants recommended for this equipment are Sentinel OPG (6700-0067-200) or Vac Kote® 37951M (0220-0091-300).**
- △ **Do not use this device in the presence of flammable anesthetics. Explosion hazard may result.**
- △ **This manual covers the operation of the venturi/ejector suction module only. Follow vacuum regulator pre-use checkout and patient set up procedures also.**
- △ **The Venturi backpack is capable of producing high levels of vacuum. Unregulated vacuum can cause injury to patients. Only use this device with an approved medical vacuum regulator.**

CAUTIONS

- △ Only trained individuals should attempt to service this device.
- △ Complete repair information is provided in the service manual.
- △ Not for field or transport use. *

* The categories of Field and Transport Use are specifically defined in ISO 10079-3 (Section 5.1.2.) "Field" means accidents or emergencies outside the hospital. "Transport" means use in ambulances, cars and airplanes. These situations may expose the equipment to uneven support, dirt, water, mechanical shock and temperature extremes. Ohio Medical suction equipment has not been tested to comply with the specific requirements of these categories.

Definitions

CAUTION: = Possibility of damage to the equipment exists.

WARNING: = Possibility of injury to the patient or the operator exists.

Note: = Provides additional information to clarify a point in the text.

Important: = Similar to a note but of greater emphasis

△ = Attention. Alerts you to a warning or caution in the text.

 = Pull knob out for I (ON)

 = Push knob in for O (OFF)

 = Supply gas inlet pressure

 = Do not use petroleum based lubricants

 = Vacuum

PTFE = Polytetrafluoroethylene

mmHg = Millimeters of Mercury

kPa = Kilopascal

mL = Milliliter

mm = Millimeter





PTS = Push-To-Set™

ISU = Intermittent Suction Unit

CVR = Continuous Vacuum Regulator

Venturi/ejector suction module identification

- ① Venturi/ejector manifold
- ② Vacuum regulator
- ③ Venturi/ejector ON/OFF switch
- ④ Supply gas probe/adaptor port
- ⑤ Exhaust
- ⑥ Positive pressure relief valve (recessed)

Equipment Set up

If the venturi/ejector suction module is purchased separately, apply a small amount of Lox-8[®] paste or Teflon[®] (PTFE) pipe sealant tape to NPT threads of the adapter on the back of the unit. Tighten adapter securely to the vacuum regulator (minimum torque of 5.4 Nm). Ensure that the unit is vertical.

- WARNINGS**
- ⚠ **Connect supply gas adapter to a 45 to 80 PSI (310 to 552 kPa) medical air or oxygen supply source only to prevent patient or operator injury and equipment damage.**
 - ⚠ **To prevent patient injury, do not block the positive pressure relief.**

- CAUTION**
- ⚠ To prevent fluid from entering the regulator and venturi/ejector suction module, use an Ohio Medical overflow safety trap and high flow bacterial suction filter. Flooding may impair operation.

High Flow Suction Filters

Hydrophilic:

Nipple

20 Pack 6730-0350-800

200 Pack 6730-0351-800


Hydrophobic:

Nipple

3 Pack 6700-0570-800 6700-0580-800

10 Pack 6700-0571-800 6700-0581-800

50 Pack 6700-0572-800 6700-0582-800


- ① Overflow safety trap
- ② High flow suction filter

Gas Consumption Rate


30 - 50 l/min.

Mode Selection

O - No vacuum is supplied to the regulator.

I - Vacuum is continuously supplied to the regulator.

Pre-Use Checkout Procedure

WARNING \triangle Perform the Pre-use Checkout Procedure before using on each patient. If there is any failure, remove from service and repair.

1. Check that the positive pressure relief valve flips open easily. You may wish to use a flat head screwdriver, wooden toothpick or silimar tool to access the relief valve flap, as long as the tool is not sharp enough to damage the flap.
2. Push the venturi ON/OFF knob to O (OFF). 
3. Turn vacuum regulator to I (ON).
Rotate the suction control knob one full turn clockwise (increasing).
Clamp the tubing to occlude the fitting port. The gauge needle should not move.
4. Turn the vacuum regulator to O (OFF)
Pull the venturi ON/OFF knob to I (ON) 
Follow the pre-use checkout procedure given in the instructions for the vacuum regulator.

Cleaning / Troubleshooting

Cleaning

Wipe all exterior surfaces with a solution of water and mild detergent after each use.

If contamination occurs, clean according to service manual recommendations.

Important: Always follow your hospital's policy for oxygen compatible cleaning; some cleaners are highly combustible in the presence of oxygen

WARNINGS  **Handle contaminated product in accordance with your hospital's infection control policy.**

 **Disinfect all suction equipment before shipment to prevent exposure to hazardous contamination.**

CAUTION  Do not steam autoclave or liquid sterilize the venturi/ejector suction module. Severe impairment to the operation of the venturi/ejector suction module will result.

Troubleshooting

If the regulator and venturi/ejector suction module do not operate, use the following:

Problem	Possible Cause	Remedy
No suction	Venturi is in the O (OFF) position	Switch to I (ON)
	Vacuum regulator is OFF 	Turn vacuum regulator on 
Reduced pressure or flow	Venturi muffler is dirty or clogged	Replace muffler
	Bacterial filter clogged	Replace bacterial filter

Important: Refer to the service manual for additional information.

This Product is sold by Ohio Medical, LLC under the warranties set forth in the following paragraphs. Such warranties are extended only with respect to the purchase of this Product directly from Ohio Medical or Ohio Medical's Authorized Dealers as new merchandise and are extended to the first Buyer thereof, other than for purpose of resale.

For a period of twelve (12) months from the date of original delivery to Buyer or to Buyer's order, but in no event for a period of more than two years from the date of original delivery by Ohio Medical, LLC to an Ohio Medical Authorized Dealer, this Product, other than its expendable parts, is warranted to be free from functional defects in materials and workmanship and to conform to the description of the Product contained in this operation manual and accompanying labels and/or inserts, provided that the same is properly operated under conditions of normal use, that regular periodic maintenance and service is performed and that replacements and repairs are made in accordance with the instructions provided. This same warranty is made for a period of thirty (30) days with respect to the expendable parts. The foregoing warranties shall not apply if the Product has been repaired other than by Ohio Medical or in accordance with written instructions provided by Ohio Medical, or altered by anyone other than Ohio Medical, or if the Product has been subject to abuse, misuse, negligence, or accident.

Ohio Medical's sole and exclusive obligation and Buyer's sole and exclusive remedy under the above warranties is limited to repairing or replacing, free of charge, at Ohio Medical's option, a Product, which is telephonically reported to the nearest Ohio Medical Regional Service Office and which, if so advised by Ohio Medical, is thereafter returned with a statement of the observed deficiency, not later than seven (7) days after the expiration date of the applicable warranty, to the designated Ohio Medical Service Office during normal business hours, transportation charges prepaid, and which, upon Ohio Medical's examination, is found not to conform with the above warranties. Ohio Medical shall not be otherwise liable for any damages including but not limited to incidental damages, consequential damages, or special damages.

There are no express or implied warranties which extend beyond the warranties hereinabove set forth. Ohio Medical, LLC makes no warranty of merchantability or fitness for a particular purpose with respect to the product or parts thereof.

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EMERGO EUROPE
Prinsessegracht 20 2514 AP The Hague The Netherlands

AUSTRALIAN SPONSOR: EMERGO AUSTRALIA
Level 20, Tower II Darling Park 201 Sussex Street Sydney, NSW 2000 Australia

© 2020 Ohio Medical, LLC.
This document contains information that is proprietary and confidential to Ohio Medical, LLC.
Use of this information is under license from Ohio Medical, LLC.
Any use other than that authorized by Ohio Medical, LLC is prohibited.

Ohio Medical and the Ohio Medical Logo are registered trademarks and Push-To-Set is a trademark of Ohio Medical, LLC.
Vac Kote is a registered trademark of the Ball Aerospace.
Lox-8 is a registered trademark of Fluoramics, Inc.
Teflon is a registered trademark of E. I. du Pont de Nemours and Company
Sentinel is a registered trademark of Sentinel Lubricants, Inc.