



VENTILATOR INTEGRATION GUIDE*

GUIDANCE ONLY - Please refer to Manufacturers Instructions For Use

Manufacturer	Model	Aerogen® integration
Acutronic	fabian HFO	Aerogen USB Controller
	fabian +nCPAP evolution	
	fabian +nCPAP evolution	
Dräger	Evita V500	Aerogen USB Controller
	Evita V300	
	Babylog VN500	
GE/Datex-Ohmeda	Engström Carestation	Standard. Dedicated cable. AG-AP1085-GE
	CARESCAPE R860	Standard. Dedicated cable. AG-AP1085-GE
Getinge (Maquet)	Servo I	Optional. Dedicated cable. AG-AP1085-MQ-I
	Servo U	Standard - Fully Integrated
	Servo N	Standard - Fully Integrated
	Servo Air	Standard - Fully Integrated
Hamilton Medical	G5	Optional. Dedicated cable. AG-AP1085-HM
	S1	Optional. Dedicated cable. AG-AP1085-HM
	C1	Aerogen USB Controller
	C2	Aerogen USB Controller
	C3	Aerogen USB Controller
	T1	Aerogen USB Controller
	MR1	Aerogen USB Controller
	C6	Optional. Dedicated cable. AG-AP1085-HM
IMT	Bellavista 1000	Aerogen USB Controller
	Bellavista neo	
	Bellavista 1000e	
	Bellavista mr	
Mindray	SV800/SV600	Aerogen USB Controller
Philips	V680	Aerogen USB Controller
ResMed	Astral	Aerogen USB Controller
SLE	SLE6000	Aerogen USB Controller